

Product Information

Version 3.1.01

A2442 DNase I (Deoxyribonuclease I)

Source: from Bovine Pancreas
Activity: >500 Kunitz units/mg
CAS #: 9003-98-9
MW: 31 kDa
Storage: Store at -20°C.

Description:

Bovine pancreatic deoxyribonuclease I (DNase I) is an endonuclease which splits phosphodiester linkages, preferentially adjacent to a pyrimidine nucleotide yielding polynucleotides with free hydroxyl group at the 3' position and phosphate group at the 5' position.

DNase I is activated by bivalent metal ions (Ca^{2+} , Mg^{2+}), and inhibited by chelating agents such EDTA and sodium dodecyl sulfate. DNase I can be stabilized against proteolytic digestion by 5mM Ca^{2+} . Optimal pH: 7-8.

Unit Definition:

One Kunitz unit will produce A260 of 0.001 per minute at pH 5.0 at 25°C using DNA, type I or III as the substrate.

Reconstitution:

Reconstitute the lyophilizate in sterile double-distilled water (10 mg/ml); further dilution with PBS (phosphate buffered saline), HBSS (Hank's Balanced Salt Solution) or medium.

Preparation Note:

1. *DNase I is sensitive to physical denaturation. Mix gently by inverting the tube. Do not vortex.*
2. For each cell type the working concentration has to be determined individually. For optimal activity the enzyme needs 5 mM Mg^{2+} .
3. Under different buffer conditions the amount of DNase required to completely digest a given amount of DNA may need to be empirically determined. For example, salt concentration >100mM will reduce DNase activity.
4. DNase I is sensitive to physical denaturation. Mix gently by inverting the tube. Do not vortex.
DNase I solution at 10 mg/ml in 150 mM NaCl may lose 10% of its activity when stored for a week at -20°C.
DNase I solution at 10 mg/ml in 150 mM NaCl may lose 20% of its activity when stored for a week at 4°C.

Enzyme Activity

Enzyme activity is assayed in the following procedure:

1) Reagent preparation

- a) 1M sodium acetate buffer with PH=5.0: Weigh 4.1g anhydrous sodium acetate, add 40ml pure water to dissolve, adjust pH=5.0 with concentrated hydrochloric acid, and finally make up to 50ml.
- b) 100 mM magnesium sulfate solution: Weigh 1.23 g of magnesium sulfate heptahydrate, dissolve it in pure water and dilute to 50 ml.
- c) 0.85% NaCl solution: Weigh 0.425g of sodium chloride, add pure water to dissolve and dilute to 50ml.
- d) DNA solution: Weigh 3.3 mg of sodium deoxyribonate type I. Add 10 ml of pure water to dissolve on ice and dissolve for at least 30 min or more until completely dissolved.
- e) Enzyme solution of DNase I (>500U/mg). Weigh 4.28mg, add 5ml 0.85% NaCl solution to dissolve.
- f) Substrate solution: 5 ml of pH=5 buffer solution, 2.5 ml of 100 ml of magnesium sulfate solution, 6 ml of DNA solution, and finally made up to 50 ml with pure water in a 50 ml volumetric flask.
- g) Blank reference solution: 5 ml of pH=5 buffer solution, 2.5 ml of 100 ml of magnesium sulfate solution, and finally made up to 50 ml with pure water in a 50 ml volumetric flask.

2) Experimental steps

Activity assay	TEST
Substrate solution	2.9
Enzyme solution	0.1

The reading is recorded every 1 minute at 260 nm, the change value of the absorbance per minute is calculated, and the maximum change value is taken to calculate the activity.

3) Calculation formula

$$CF = R_{std} / T_{std}$$

$$Activity = K * 3 * 5 / 0.001 * v * m$$

$$Corrected\ activity = activity * CF$$

Where: CF- activity correction factor

R_{std} - standard nominal activity, U/mg

T_{std} - standard test activity, U/mg

V - the volume of the enzyme solution added, ml

3- Total volume of reaction solution, ml

5-Preparation of the volume of the enzyme solution, ml

0.001-A260nm per 0.001 absorbance is defined as one unit activity

This material is for laboratory research purpose and/or in vitro use only and is not to be used in humans or animals.