

Product Information

Version 1.1.07

A3806– RNase A, Pancreas, >70U/MG

Alternative names: Ribonuclease I; RNase I; Pancreatic ribonuclease; RNase A; Endoribonuclease I

Systematic name: Ribonuclease 3'-pyrimidino-oligonucleotidohydrolase

Source: Bovine pancreas. The raw material used is of Australian and New Zealand origin, collected from abattoirs which are under veterinary control. Furthermore, the product was exposed to a pH of 4.5 for 2 days and a temperature of 80 °C for approx. 5 minutes during purification.

Form: An essentially protease free, chromatographically prepared freeze dried material

Unit Definition: One unit causing the hydrolysis of RNA at a rate such that k (velocity constant) equals unity (Kunitz units) at 25 °C and pH 5.0.

Activity >70 Kunitz/mg material (RNase A content approximately 70%)

Solubility: Dissolves readily at 2 mg/ml in analytical grade water to give a clear colorless solution.

Product Description

A major application for Ribonuclease A (RNase A) is the removal of RNA from preparations of plasmid DNA. In this application, the presence of DNase activity as an impurity is a concern. The boiling-water bath method used to eliminate contaminating DNase activity has proven unreliable. RNase A is an endoribonuclease that attacks at the 3' phosphate of a pyrimidine nucleotide. The sequence of pG-pG-pC-pA-pG will be cleaved to give pG-pG-pCp and A-pG. The highest activity is exhibited with single stranded RNA. RNase A is a single chain polypeptide containing 4 disulfide bridges. In contrast to RNase B, it is not a glycoprotein. RNase A can be inhibited by alkylation of His12 or His119, which are present in the active site of the enzyme. Activators of RNase A include potassium and sodium salts.

Molecular weight: 13.7 kDa (amino acid sequence)

Extinction coefficient: E1% = 7.0 (280 nm)

Isoelectric point: pI = 9.6

Optimal temperature: 60 °C (activity range of 15-70 °C)

Optimal pH: 7.6 (activity range of 6-10)

Inhibitors: ribonuclease inhibitor

Note: RNase A is stable to both heat and detergents. In addition, it adsorbs strongly to glass. Scrupulous precautions are necessary to ensure that residues of RNase A do not cause artifacts in processes requiring intact RNA.

Application

- RNase protection assays
- Remove unspecifically bound RNA
- Analysis of RNA sequences
- Hydrolyze RNA contained in protein samples
- Purification of DNA

Pancreatic RNase A specifically cleaves at the 3'-side of pyrimidine (uracil or cytosine) phosphate bonds. Ribonucleases do not hydrolyze DNA, because the DNA lacks 2'-OH groups essential for the formation of cyclic intermediates. RNase can also hydrolyze RNA from protein samples.

Storage Condition

-20 °C. Protect from moisture. Allow to come to room temperature before opening.

RNase A is a very stable enzyme and solutions have been reported to withstand temperatures up to 100 °C. At 100 °C, an RNase A solution is most stable between pH 2.0 and 4.5.

Preparation Instructions

A stock solution of 10 mg/ml RNase A could be prepared by dissolving 100mg of RNase A in 9.925ml of 0.01M NaOAc (pH 5.2), plus 0.075ml of 1M Tris (not pH adjusted). Please verify final pH to be neutral.

References

1. Smyth, D.G. et. al., J. Biol. Chem., 238, 227-234 (1963).
2. Pace, C.N. et. al., Protein Science, 4, 2411-2423 (1995).
3. Tanford, C., and Hauenstein, J. D., J. Am. Chem. Soc., 78, 5287-5291 (1956).
4. Plummer, T.H., and Hirs, C.H.W., J. Biol. Chem., 238, 1396-1397 (1963).
5. Heinrikson, R.L. et. al., J. Biol. Chem., 240, 2921- 2934 (1965).
6. Schomberg, D., and Salzmann, M., Enzyme Handbook, Vol. 3, 1-3 under E.C. 3.1.27.5 (1990).
7. Burrell, M.M., Enzymes of Molecular Biology, Vol. 16, 263-270 (1993).
8. Crestfield, A.M. et. al., J. Biol. Chem., 238, 618-621 (1963).
9. Sambrook, J. et. al., Molecular Cloning, A Laboratory Manual, 2nd ed., 1.51 (1989).
10. Sambrook, J. et. al., Molecular Cloning, A Laboratory Manual, 2nd ed., B.17.

**** This product is for laboratory research purpose only. Not for human or animal diagnostic and therapeutic use.***