



Tel: (519) 489-7195, (800) 836-8089 Fax: (519) 231-0140, (877) 221-3515

Email: info@biomatik.com https://www.biomatik.com

User Manual

Catalog Number: EKC42224

Product Name: Human carbonhydrate antigen 19-9 (CA19-9) ELISA kit

Detection Range:10 U/mL-150 U/mL

Intended Use: For quantitative determination of Human carbonhydrate antigen

19-9 (CA19-9) concentrations in serum, plasma.

Precautions: For research use only. Not for use in diagnostic procedures.

Manual Version: 202301V1

Storage:

| Unopened kit | 6 months when stored at 2 - 8°C. |
|--------------|---|
| Opened Kit | May be stored up to 1 month at 2 - 8°C. Keep it in sealed aluminum foil bag and avoid moisture. |

The product manual may be updated as a result of continuous improvements.

Always refer to the hard copy manual included in the kit for your experiment.





Tel: (519) 489-7195, (800) 836-8089 Fax: (519) 231-0140, (877) 221-3515

Email: info@biomatik.com https://www.biomatik.com

Kit Components

| Reagents | Quantity | | |
|--|-------------|--|--|
| Assay plate (12 x 8 coated Microwells) | 1(96 wells) | | |
| Standard | 5 x 0.5 ml | | |
| HRP-conjugate | 1 x 6 ml | | |
| Antibody | 1 x 6 ml | | |
| Wash Buffer (20 x concentrate) | 1 x 15 ml | | |
| Substrate A | 1 x 7 ml | | |
| Substrate B | 1 x 7 ml | | |
| Stop Solution | 1 x 7 ml | | |
| Adhesive Strip (For 96 wells) | 4 | | |
| Instruction manual | 1 | | |

Standard Concentration

| Standard | S0 | S1 | S2 | S3 | S4 |
|-------------------------|----|----|----|-----|-----|
| Concentration (U/ml) | 0 | 10 | 40 | 100 | 150 |

Working Principle

This assay employs the quantitative sandwich enzyme immunoassay technique. Antibody specific for CA19-9 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells with a Horseradish Peroxidase (HRP) conjugated antibody specific for CA19-9. Following a wash to remove any unbound reagent, a substrate solution is added to the wells and color develops in proportion to the amount of CA19-9 bound in the initial step. Measure the intensity of the color when the color stopped developing.

Tel: (519) 489-7195, (800) 836-8089 Fax: (519) 231-0140, (877) 221-3515

Email: info@biomatik.com

https://www.biomatik.com

Sensitivity

BIOMATIK

The minimum detectable dose of human CA19-9 is less than 3.8 U/mL.

The sensitivity of this assay, or Lower Limit of Detection (LLD) was

defined as the lowest human CA19-9 concentration that could be

differentiated from zero. It was determined by the mean O.D value of 20

replicates of the zero standard added by their three standard deviations.

Specificity

This assay has high sensitivity and excellent specificity. No significant

cross-reactivity or interference between human CA19-9 and analogues

was observed. Limited by current skills and knowledge, it is impossible

for us to complete the cross-reactivity detection between human CA19-9

and all the analogues, therefore, cross reaction may still exist.

Precision

Intra-assay Precision (Precision within an assay): CV%<15%

Three samples of known concentration were tested twenty times on one

plate to assess.

Inter-assay Precision (Precision between assays): CV%<20%

Three samples of known concentration were tested in twenty assays to

assess.

Tel: (519) 489-7195, (800) 836-8089 Fax: (519) 231-0140, (877) 221-3515

Email: info@biomatik.com

https://www.biomatik.com

Other Supplies Required

BIOMATIK

Microplate reader capable of measuring absorbance at 450 nm,

with the correction wavelength set at 540 nm or 570 nm.

An incubator which can provide stable incubation conditions up to

37°C±0.5°C.

Squirt bottle, manifold dispenser, or automated microplate washer.

Absorbent paper for blotting the microtiter plate.

100ml and 500ml graduated cylinders.

Deionized or distilled water.

Pipettes and pipette tips.

Test tubes for dilution.

Sample Collection & Storage

Serum Use a serum separator tube (SST) and allow samples to

clot for two hours at room temperature or overnight at 4°C before

centrifugation at 1000 ×g for 15 mins. Remove serum and assay

immediately or aliquot and store samples at -20°C or -80°C. Avoid

repeated freeze-thaw cycles.

Tel: (519) 489-7195, (800) 836-8089 Fax: (519) 231-0140, (877) 221-3515

Email: info@biomatik.com

https://www.biomatik.com

Plasma Collect plasma using EDTA, or heparin as

anticoagulant. Centrifuge for 15 minutes at 1000 ×g at 2-8°C within 30

minutes of collection. Assav immediately or aliquot and store samples at -

20°C or -80°C. Avoid repeated freeze-thaw cycles.

BIOMATIK

Notes:

1. Biomatik is only responsible for the kit itself, not for the samples

consumed during the assay. The user need to calculate the possible

amount of the samples to be used in the whole test. Please reserve

sufficient samples in advance.

2. Samples to be used within 5 days may be stored at 2-8°C, otherwise,

samples must be stored at -20°C (≤1month) or -80°C (≤2month) to avoid

contamination and loss of bioactivity.

3. Grossly hemolyzed samples are not suitable to use.

4. It would be necessary to run a preliminary experiment for validation, if

the samples are not indicated in the manual.

5. Please predict the concentration before assaying. If results were not

within the range of the standard curve, users would need to estimate the

optimal sample dilutions for their particular experiments.

Tel: (519) 489-7195, (800) 836-8089

Fax: (519) 231-0140, (877) 221-3515

Email: info@biomatik.com https://www.biomatik.com

BIOMATIK

6. Tissue or cell extraction samples prepared by chemical lysis buffer

may cause unexpected ELISA results due to the impacts of certain

chemicals

7. Considering the possibility of mismatch between antigen from other

resources and antibody in our kits (e.g., antibody targets conformational

epitope rather than linear epitope), some native or recombinant proteins

from other manufacturers may not be recognized by our products.

8. Due to factors including cell viability, cell number and sampling time.

samples from cell culture supernatant may not be detected by the kit.

9. Recommend to use fresh samples for the test. If you store samples for

long time, protein degradation and denaturalization may occur in those

samples and finally lead to wrong results.

Reagents Preparation

Kindly use graduated containers to prepare the reagent.

Please don't prepare the reagent directly in the Diluent vials in the

kit.

Bring all reagents to room temperature (18-25°C) before use for 30

mins

Distilled Water is recommended. Contaminated water or container for

reagents preparation will affect the test result.

Tel: (519) 489-7195, (800) 836-8089

Fax: (519) 231-0140, (877) 221-3515

Email: info@biomatik.com https://www.biomatik.com

Wash Buffer(1x)- If crystals have formed in the concentrate, warm up to

room temperature and mix gently until the crystals have completely

dissolved. Dilute 15 ml of Wash Buffer Concentrate (20 x) into deionized

or distilled water to prepare 300 ml of Wash Buffer (1 x).

HRP-conjugate mixed solution - V_{Antibody} -V_{HRP-conjugate} =1:10 (eg: 0.1ml

Antibody + 1ml HRP-conjugate, shake well.

BIOMATIK

Key Notes

Do not mix or substitute reagents with those from other lots or

sources.

If samples generate higher values than the highest standard, dilute

the samples with Sample Diluent and repeat the assay.

Any variation in Sample Diluent, operator, pipetting technique,

washing technique, incubation time or temperature, and kit age can

cause variation in binding.

This assay is designed to eliminate interference by soluble

receptors, binding proteins, and other factors present in biological

samples. Until all factors have been tested in the Immunoassay, the

possibility of interference cannot be excluded.

Tel: (519) 489-7195, (800) 836-8089 Fax: (519) 231-0140, (877) 221-3515

Email: info@biomatik.com https://www.biomatik.com

Precautions

The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face,

and clothing protection when using this material.

BIOMATIK

Assay Procedures

Bring all reagents and samples to room temperature before use.

Centrifuge the sample again after thawing before the assay. It is

recommended to assay all samples and standards in duplicate.

1. Prepare all reagents, working standards and samples as directed in

the previous sections.

2. Determine the number of wells to be used and put any remaining wells

and the desiccant back into the pouch and seal the ziploc, store unused

wells at 4°C.

3. Set a Blank well without any solution.

4. Add 50µl Standard or Sample per well.

5. Add 50µl **HRP-conjugate** to each well(Not to Blank well). Mix well and

then incubate at 37°C for 1 hour.

6. Aspirate each well and wash, repeating the process two times for a

total of three washes. Wash by filling each well with Wash Buffer (200µl)

using a squirt bottle, multi-channel pipette, manifold dispenser, or

autowasher, and let it stand for 10 seconds, complete removal of liquid at



Tel: (519) 489-7195, (800) 836-8089 Fax: (519) 231-0140, (877) 221-3515

Email: info@biomatik.com

https://www.biomatik.com

BIOMATIK

each step is essential to good performance. After the last wash, remove

any remaining Wash Buffer by aspirating or decanting. Invert the plate

and blot it against clean paper towels.

7. Add 50µl Substrate A and 50µl Substrate B to each well, mix well.

Incubate at 37°C for 15 mins. Keeping the plate away from drafts and

other temperature fluctuations in the dark.

8. Add 50µl Stop Solution to each well, gently tap the plate to ensure

thorough mixing.

9. Determine the optical density of each well within 10 mins, using a

microplate reader set to 450 nm.

Notes:

1 The final experimental results will be closely related to validity of

the products, operation skills of the end users and the experimental

environments.

2. Samples or reagents addition: Please use the freshly prepared

Standard. Please carefully add samples to wells and mix gently to avoid

foaming. Do not touch the well wall as possible. For each step in the

procedure, total dispensing time for addition of reagents or samples to

the assay plate should not exceed 10 mins. This will ensure equal

elapsed time for each pipetting step, without interruption. Duplication of



BIOMATIK

Tel: (519) 489-7195, (800) 836-8089

Fax: (519) 231-0140, (877) 221-3515 Email: info@biomatik.com https://www.biomatik.com

all standards and specimens, although not required, is recommended. To

avoid cross-contamination, change pipette tips between additions of each

standard level, between sample additions, and between reagent additions.

Also, use separate reservoirs for each reagent.

3. Incubation: To ensure accurate results, proper adhesion of plate

sealers during incubation steps is necessary. Do not allow wells to sit

uncovered for extended periods between incubation steps. Once

reagents have been added to the well strips. DO NOT let the strips DRY

at any time during the assay. Incubation time and temperature must be

observed.

Washing: The wash procedure is critical. Complete removal of 4.

liquid at each step is essential to good performance. After the last wash,

remove any remaining Wash Solution by aspirating or decanting and

remove any drop of water and fingerprint on the bottom of the plate.

Insufficient washing will result in poor precision and falsely elevated

absorbance reading. When using an automated plate washer, adding a

30 second soak period following the addition of wash buffer, and/or

rotating the plate 180 degrees between wash steps may improve assay

precision.





Tel: (519) 489-7195, (800) 836-8089 Fax: (519) 231-0140, (877) 221-3515

Email: info@biomatik.com

https://www.biomatik.com

5. Controlling of reaction time: Observe the change of color after adding TMB Substrate (e.g. observation once every 10 mins), TMB Substrate should change from colorless or light blue to gradations of blue. If the color is too deep, add Stop Solution in advance to avoid excessively strong reaction which will result in inaccurate absorbance reading.

- 6. TMB Substrate is easily contaminated. TMB Substrate should remain colorless or light blue until added to the plate. Please protect it from light.
- Stop Solution should be added to the plate in the same order as 7. the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that Stop Solution has not mixed thoroughly with TMB Substrate.

Tel: (519) 489-7195, (800) 836-8089 Fax: (519) 231-0140. (877) 221-3515

Email: info@biomatik.com

https://www.biomatik.com

Calculation of Results

BIOMATIK

Using the professional soft "Curve Expert" to make a standard

curve is recommended, which can be downloaded from our web.

Average the duplicate readings for each standard and sample and

subtract the average optical density of Blank.

Create a standard curve by reducing the data using computer software

capable of generating a four parameter logistic (4-PL) curve-fit. As an

alternative, construct a standard curve by plotting the mean absorbance

for each standard on the x-axis against the concentration on the y-axis

and draw a best fit curve through the points on the graph. The data may

be linearized by plotting the log of the CA19-9 concentrations versus the

log of the O.D. and the best fit line can be determined by regression

analysis. This procedure will produce an adequate but less precise fit of

the data.

If samples have been diluted, the concentration read from the standard

curve must be multiplied by the dilution factor.