

Certificate of Quality

Product Name: dNTP Mix, 25mM Each

A solution contains 25mM dATP, 25mM dCTP, 25mM dGTP, and 25mM dTTP.

Product Code:	A2447	CAS #:	
Grade:	Biotechnology	Formula:	
Appearance:	Clear, colorless liquid	MW:	
Reassay Date:	2 years after shipping date	Lot#:	See product label
Storage:	-20°C. Allow to come to room temperature before opening. avoid freeze-thaw cycles.		

Quality Assurance Statement: Functionally tested in PCR with Taq DNA polymerase, Pfu DNA polymerase and T7 DNA polymerase.

TEST	SPECIFICATION	RESULT
DNase & RNase	None	None Detected
pH @25°C	7.0	Passes

Data Verified by M. He. Verification Date: 2020-01-06

Limited Usage and Warranty: For Research Use Only. Not for use in diagnostic procedures. Not for human or animal drug or food use.

Biomatik warrants materials of said quality at the time of sale. It is the sole responsibility of the customer to determine the adequacy of all materials for any intended or specific purpose of use. Biomatik's sole obligation is to replace the material up to the extent of the purchase price. This warranty applies only to products in original packaging and does not apply to a product which has been tampered with or altered in any way or which has been misused or damaged by accident or negligence. No other warranties of any kind, expressed or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Biomatik. Biomatik shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.