

Certificate of Quality

Product Name: Trypsin Inhibitor, Soybean-Purified

A serine protease inhibitor specific for trypsin and trypsin-like proteases.

Product Code: A5107 CAS #: 9035-81-8
Grade: High Purity Formula: N/A
Appearance: Off-white to beige lyophilized powder MW: N/A
Reassay Date: 2 years after shipping date Lot#: See product label
Storage: -20°C. Protect from moisture. Allow to come to room temperature before opening, avoid freeze-thaw cycles.

TEST	SPECIFICATION	RESULT
Activity (material)	≥ 10,000 BAEE units of inhibition/mg	Passes
Solubility	Passes	Passes
Protein Conc	0.872 mg protein/material	Passes

Data Verified by M. He. Verification Date: 2020-01-06

Limited Usage and Warranty: For Research Use Only. Not for use in diagnostic procedures. Not for human or animal drug or food use. Biomatik warrants materials of said quality at the time of sale. It is the sole responsibility of the customer to determine the adequacy of all materials for any intended or specific purpose of use. Biomatik's sole obligation is to replace the material up to the extent of the purchase price. This warranty applies only to products in original packaging and does not apply to a product which has been tampered with or altered in any way or which has been misused or damaged by accident or negligence. No other warranties of any kind, expressed or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Biomatik. Biomatik shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.