

Certificate of Quality

Product Name: 1Kb DNA Ladder

Product Code:	M7508	CAS #:	N/A
Grade:	Biotechnology	Formula:	N/A
Appearance:		MW:	N/A
Reassay Date:	2 years after shipping date	Lot#:	See product label
Storage:	-20 °C in a non-frost-free freezer		

Quality Assurance Statement: 1Kb DNA ladder (Quantitative) is ideal for determining the size of double-stranded DNA from 500 to 1000 base pairs. The ladder consists of 10 linear double-stranded DNA fragments, chromatography-purified. The intensity of the 2000bp and 5000bp has been increased to serve as a reference for easy identification. All fragments are precisely quantified and mixed during the production. For 5µl loading, all fragments except 2000bp and 5000bp are at 40ng. The 2000bp and 5000bp fragments are at 100ng. This ladder is pre-mixed with loading dye and is ready to use.

TEST	SPECIFICATION	RESULT
------	---------------	--------

Data Verified by M. He. Verification Date: 2020-01-06

Limited Usage and Warranty: For Research Use Only. Not for use in diagnostic procedures. Not for human or animal drug or food use. Biomatik warrants materials of said quality at the time of sale. It is the sole responsibility of the customer to determine the adequacy of all materials for any intended or specific purpose of use. Biomatik's sole obligation is to replace the material up to the extent of the purchase price. This warranty applies only to products in original packaging and does not apply to a product which has been tampered with or altered in any way or which has been misused or damaged by accident or negligence. No other warranties of any kind, expressed or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Biomatik. Biomatik shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.