



Biomatik
Tel: (519) 489-7195, (800) 836-8089
Fax: (519) 231-0140, (877) 221-3515
Email: info@biomatik.com
http://www.biomatik.com

Product Information

Version 1.1.03

A3614– pH Standard Kit

Kit Contents	Catalog #	Grade	Size
pH Standard, 10.0	A3611	Laboratory	500ML
pH Standard 4.0	A3612	Laboratory	500ML
pH Standard 7.0	A3613	Laboratory	500ML

Description

Kit contains 500ml each of pH 4.0, pH 7.0 and pH 10.0 standards, Color coded, accuracy to NIST standards.

pH Standard, 10.0

Catalog Number: A3611	Formula: N/A
CAS Number: N/A	Molecular Weight: N/A
Grade: Laboratory	Appearance: Blue, haze-free solution

Specifications: pH @25C: 9.97 – 10.03

pH Standard, 4.0

Catalog Number: A3612	Formula: N/A
CAS Number: N/A	Molecular Weight: N/A
Grade: Laboratory	Appearance: Pink, haze-free solution

Specifications: pH @25C: 3.97 - 4.03

pH Standard, 7.0

Catalog Number: A3613	Formula: N/A
CAS Number: N/A	Molecular Weight: N/A
Grade: Laboratory	Appearance: Yellow, haze-free solution

Specifications: pH @25C: 6.97 - 7.03

Storage Condition

Stable at Room temperature.

**** This product is for laboratory research purpose only. Not for human or animal diagnostic and therapeutic use.***