

## Anti-SARS-CoV-2 antibody

**Catalog No.:** CAD10001

**Application:** SARS-CoV-2 Kit

**Description:** Cultured in vitro, without contamination of mouse antibodies

**Source:** Mammalian cell

**Species reactivity:** React with S protein

**Isotype:** Monoclonal IgG

**Form:** Clear and colorless liquid without flocculent precipitate

**Purity:** > 98%, purified by Protein A/G chromatography

**Buffer:** 1 × PBS, pH7.4

**Storage:** Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Store at 2 - 8°C for up to 12 months or at -20 °C for up to 24 months.