

anti- AKAP3 antibody

Product Information

Catalog No.:	CAF50015
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May function as a regulator of both motility-and head-associated functions such as capacitation and the acrosome reaction.

Immunogen information


Immunogen:	A kinase(PRKA) anchor protein 3
Synonyms:	A kinase anchor protein 3, AKAP 110, AKAP 3, AKAP110, AKAP3, Cancer/testis antigen 82, CT82, Fibrousheathin 1, Fibrousheathin I, FSP95, PRKA3, SOB1, Sperm oocyte binding protein
Observed MW:	95 kDa
Uniprot ID :	O75969

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:10000; IHC: 1:20-200
Image:	

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Immunohistochemistry of paraffin-embedded human testis using CAF50015(AKAP3 antibody) at dilution of 1:100



mouse testis tissue were subjected to SDS PAGE followed by western blot with CAF50015(AKAP3 antibody) at dilution of 1:5000