

anti- Caspase 2 antibody

Product Information

Catalog No.:	CAF50060
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information


Immunogen:	caspase 2, apoptosis-related cysteine peptidase
Synonyms:	CASP 2, CASP2, Caspase 2, ICH 1L, ICH 1L/1S, ICH1, NEDD 2, NEDD2, Protease ICH 1
Observed MW:	50 kDa
Uniprot ID :	P42575

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Immunohistochemistry of paraffin-embedded rat kidney tissue slide using CAF50060(Caspase 2 Antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with CAF50060(CASP2 antibody) at dilution of 1:1000