

anti- Cyclin A2 antibody

Product Information


Catalog No.:	CAF50080
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20 $^{\circ}$ C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	cyclin A2
Synonyms:	CCN1, CCNA, CCNA2, Cyclin A, Cyclin A2
Observed MW:	48-52 kDa
UniprotID :	P20248

Application

Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB, IF, IHC		
Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:50-1:200			
Image:			



Immunohistochemistry of paraffin-embedded human squamous skin carcinoma using CAF50080(cyclin A2 antibody) at dilution of 1:100



Biomatik Tel: (519) 489-7195, (800) 836-8089 Fax: (519) 231-0140, (877) 221-3515 Email: info@ biomatik.com https://www.biomatik.com

