

## anti- CYP7B1 antibody

### Product Information

Catalog No.:	CAF50083
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	cytochrome P450, family 7, subfamily B, polypeptide 1
Synonyms:	CBAS3, CP7B, CYP7B1, Cytochrome P450 7B1, Oxysterol 7 alpha hydroxylase, SPG5A
Observed MW:	50-55 kDa
Uniprot ID :	O75881

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200; IF: 1:20-1:200
Image:	

**This Antibody is for Research Use Only. Not for Diagnostic Procedures.**

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using CAF50083(CYP7B1 Antibody) at dilution of 1:50



A431 cells were subjected to SDS PAGE followed by western blot with CAF50083(CYP7B1 Antibody) at dilution of 1:1000

SAMPLE

ONLY

**This Antibody is for Research Use Only. Not for Diagnostic Procedures.**

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.