

anti- G-CSF antibody

Product Information

Catalog No.:	CAF50129
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	colony stimulating factor 3(granulocyte)
Synonyms:	C17orf33, CSF3, Filgrastim, G CSF, GCSF, G-CSF, Lenograstim, Pluripoietin
Observed MW:	22 kDa
Uniprot ID :	P09919

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using CAF50129(G-CSF Antibody) at dilution of 1:50



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with CAF50129(G-CSF Antibody) at dilution of 1:1000

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.