

anti- MAGEA10 antibody

Product Information

Catalog No.:	CAF50209
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Not known, though may play a role in embryonal development and tumor transformation or aspects of tumor progression.

Immunogen information


Immunogen:	melanoma antigen family A, 10
Synonyms:	MAGE10
Observed MW:	41kd
Uniprot ID :	P43363

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Jurkat cells were subjected to SDS PAGE followed by western blot with CAF50209(MAGEA10 antibody) at dilution of 1:500

ONLY

SAMPLE

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.