

## anti- MAGEA4 antibody

### Product Information

Catalog No.:	CAF50211
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Not known, though may play a role in embryonal development and tumor transformation or aspects of tumor progression.

### Immunogen information


Immunogen:	melanoma antigen family A, 4
Synonyms:	MAGE4
Observed MW:	34 kDa
Uniprot ID :	P43358

### Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

**This Antibody is for Research Use Only. Not for Diagnostic Procedures.**

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Immunohistochemistry of paraffin-embedded human breast cancer using CAF50211(MAGEA4 antibody) at dilution of 1:50



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with CAF50211(MAGEA4 Antibody) at dilution of 1:500