

anti- MYL4 antibody

Product Information

Catalog No.:	CAF50231
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Regulatory light chain of myosin. Does not bind calcium.

Immunogen information

Immunogen:	myosin, light chain 4, alkali; atrial, embryonic
Synonyms:	ALC1, AMLC, GT1, MLC1, MYL4, Myosin light chain 4, PRO1957
Observed MW:	22kd
Uniprot ID :	P12829

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



human liver tissue were subjected to SDS PAGE followed by western blot with CAF50231(MYL4 antibody) at dilution of 1:1000

ONLY

SAMPLE

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.