

BiomatikTel: (519) 489-7195, (800) 836-8089
Fax: (519) 231-0140, (877) 221-3515
Email: info@ biomatik.com https://www.biomatik.com

anti- NXF1 antibody

Product Information

Catalog No.: CAF50240

Size: 100µg Form: liquid

Purification: Immunogen affinity purified

≥95% as determined by SDS-PAGE Purity:

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 Storage:

months (Avoid repeated freeze / thaw cycles.)

Background

Involved in the nuclear export of mRNA species bearing retroviral constitutive transport elements(CTE) and in the export of mRNA from the nucleus to the cytoplasm(TAP/NFX1 pathway). The NXF1-NXT1 heterodimer is involved in the export of HSP70 mRNA in conjunction with ALYREF/THOC4 and THOC5 components of the TREX complex. ALYREF/THOC4-bound mRNA is thought to be transferred to the NXF1-NXT1 heterodimer for export.

Immunogen information

nuclear RNA export factor 1 Immunogen:

MEX67, mRNA export factor TAP, nuclear RNA export factor 1, NXF1, Synonyms:

TAP, Tip associated protein, Tip associating protein

Observed MW: 66-70 kDa Q9UBU9 Uniprot ID:

Application

Reactivity: Human, Mouse, Rat


Tested Application: ELISA, WB, IF


Recommended dilution: WB: 1:500-1:1000; IF: 1:20-1:200


Image:



Biomatik Tel: (519) 489-7195, (800) 836-8089 Fax: (519) 231-0140, (877) 221-3515 Email: info@ biomatik.com https://www.biomatik.com







Hela cells were subjected to SDS PAGE followed by western blot with CAF50240(NXF1 Antibody) at dilution of 1:600