

anti- RPP21 antibody

Product Information

Catalog No.:	CAF50303
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20 $^{\circ}$ C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Component of ribonuclease P, a protein complex that generates mature tRNA molecules by cleaving their 5'-ends.

Immunogen information

Immunogen:	ribonuclease P/MRP 21kDa subunit
Synonyms:	C6orf135, CAT60, Ribonucleoprotein V, RNaseP protein p21, RPP21
Observed MW:	18 kDa
Uniprot ID :	Q9H633

Application

Reactivity:	Human, Mouse, Rat			
Tested Application:	ELISA, WB			
Recommended dilution:WB: 1:500-1:1000				
Image:				



Biomatik Tel: (519) 489-7195, (800) 836-8089 Fax: (519) 231-0140, (877) 221-3515 Email: info@ biomatik.com https://www.biomatik.com

