

## anti- SKIL antibody

### Product Information

Catalog No.:	CAF50329
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May have regulatory role in cell division or differentiation in response to extracellular signals.

### Immunogen information

Immunogen:	SKI-like oncogene
Synonyms:	SNO
Observed MW:	77 kDa
Uniprot ID :	P12757



### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	

**This Antibody is for Research Use Only. Not for Diagnostic Procedures.**

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.

IP Result of anti-SKIL (IP:CAF50329, 4ug;  
Detection:CAF50329 1:500) with HepG2 cells  
lysate 1600ug.



A431 cells were subjected to SDS PAGE followed  
by western blot with CAF50329(SKIL antibody) at  
dilution of 1:300

**This Antibody is for Research Use Only. Not for Diagnostic Procedures.**

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.