

# anti- STAT6 antibody

### **Product Information**

Catalog No.:	CAF50343
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20 $^{\circ}$ C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Carries out a dual function: signal transduction and activation of transcription. Involved in IL4/interleukin-4-and IL3/interleukin-3-mediated signaling.

#### **Immunogen information**

Immunogen:	signal transducer and activator of transcription 6, interleukin-4 induced
Synonyms:	D12S1644, IL 4 STAT, STAT6, STAT6B, STAT6C
Observed MW:	100-110 kDa
Uniprot ID :	P42226

## Application

Reactivity:	Human, Mouse, Rat			
Tested Application:	ELISA, WB, IHC			
Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200				
Image:				



Biomatik Tel: (519) 489-7195, (800) 836-8089 Fax: (519) 231-0140, (877) 221-3515 Email: info@ biomatik.com https://www.biomatik.com

