

## anti- TACC1 antibody

### Product Information

Catalog No.:	CAF50349
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Likely involved in the processes that promote cell division prior to the formation of differentiated tissues.

### Immunogen information


Immunogen:	transforming, acidic coiled-coil containing protein 1
Synonyms:	KIAA1103
Observed MW:	110-130 kDa
Uniprot ID :	O75410

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	

**This Antibody is for Research Use Only. Not for Diagnostic Procedures.**

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



MCF7 cells were subjected to SDS PAGE followed by western blot with CAF50349(TACC1 antibody) at dilution of 1:1000

ONLY

SAMPLE

**This Antibody is for Research Use Only. Not for Diagnostic Procedures.**

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.