

anti- UQCR antibody

Product Information

Catalog No.:	CAF50368
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information


Immunogen:	ubiquinol-cytochrome c reductase, 6.4kDa subunit
Synonyms:	0710008D09Rik, Complex III subunit 10, Complex III subunit XI, Cytochrome b c1 complex subunit 10, QCR10, UQCR, UQCR11
Observed MW:	6 kDa
Uniprot ID :	O14957

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:20-1:200
Image:	

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Immunohistochemical of paraffin-embedded human heart using CAF50368(UQCR antibody) at dilution of 1:50

ONLY

SAMPLE

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.