

anti- WFS1 antibody

Product Information

Catalog No.:	CAF50374
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Participates in the regulation of cellular Ca(2+) homeostasis, at least partly, by modulating the filling state of the endoplasmic reticulum Ca(2+) store.

Immunogen information


Immunogen:	Wolfram syndrome 1(wolframin)
Synonyms:	FLJ51211, WFRS, WFS, WFS1, Wolfram syndrome 1(wolframin), WOLFRAMIN
Observed MW:	100 kDa
Uniprot ID :	O76024

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, IF, WB
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:500; IF: 1:50-1:200
Image:	

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using CAF50374(WFS1 Antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with CAF50374(WFS1 Antibody) at dilution of 1:600

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.