

## anti- JAK1 antibody

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | CAF50395  |
| Size:         | 100µg   |
| Form:         | liquid  |
| Purification: | Immunogen affinity purified   |
| Purity:       | ≥95% as determined by SDS-PAGE  |
| Host:         | Rabbit  |
| Clonality:    | polyclonal  |
| Clone ID:     | None  |
| IsoType:      | IgG   |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Immunogen information


|              |   |
|--------------|---|
| Immunogen:   | Janus kinase 1  |
| Synonyms:    | JAK 1, JAK1, JAK1A, JAK1B, Janus kinase 1, JTK3, Tyrosine protein kinase JAK1 |
| Observed MW: | None  |
| Uniprot ID : | P23458  |

### Application

|                            |              |
|----------------------------|--------------|
| Reactivity:                | Human        |
| Tested Application:        | ELISA, IHC   |
| Recommended dilution: IHC: | 1:50 - 1:200 |
| Image:                     |              |

**This Antibody is for Research Use Only. Not for Diagnostic Procedures.**

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Immunohistochemistry of paraffin-embedded human intestinal cancer tissue slide using CAF50395( JAK1 antibody) at dilution of 1:100

ONLY

SAMPLE

**This Antibody is for Research Use Only. Not for Diagnostic Procedures.**

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.