

anti- CDR2 antibody

Product Information

Catalog No.:	CAF50397
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	cerebellar degeneration-related protein 2
Synonyms:	CDR2, CDR62, PCD17, Yo
Observed MW:	None
Uniprot ID :	Q01850

Application

Reactivity:	Human, rat
Tested Application:	IHC
Recommended dilution:	IHC: 1:50 - 1:200
Image:	

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Immunohistochemistry of paraffin-embedded rat brain tissue slide using CAF50397(CDR2 antibody) at dilution of 1:100

ONLY

SAMPLE

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.