

anti- ACTN3 antibody

Product Information

Catalog No.:	CAF50440
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.

Immunogen information


Immunogen:	actinin, alpha 3
Synonyms:	Actinin, alpha 3, ACTN3, Alpha actinin 3, F actin cross linking protein
Observed MW:	103 kDa
Uniprot ID :	Q08043

Application


Reactivity:	Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Immunohistochemistry of paraffin-embedded rat heart tissue slide using CAF50440(ACTN3 Antibody) at dilution of 1:200



rat skeletal muscle tissue were subjected to SDS PAGE followed by western blot with CAF50440(ACTN3 Antibody) at dilution of 1:1000

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.