

## anti- PE2R4 antibody

### Product Information

Catalog No.:	CAF50450
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	PGE2 receptor EP4 subtype
Synonyms:	PTGER4, EP4, EP4R, prostaglandin E receptor 4, PE2R4
Observed MW:	50 kDa, 22 kDa
Uniprot ID :	P35408

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	

**This Antibody is for Research Use Only. Not for Diagnostic Procedures.**

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Raji cells were subjected to SDS PAGE followed by western blot with CAF50450( PE2R4 Antibody) at dilution of 1:1000

ONLY

SAMPLE

**This Antibody is for Research Use Only. Not for Diagnostic Procedures.**

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.