

anti- ABC3A antibody

Product Information

Catalog No.:	CAF50452
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	APOBEC-3A
Synonyms:	APOBEC3A, A3A, ARP3, PHRBN, bK150C2.1, DNA dC->dU-editing enzyme APOBEC-3A
Observed MW:	26 kDa
Uniprot ID :	P31941

Application


Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:500
Image:	

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using CAF50452(ABC3A Antibody) at dilution of 1:200



rat testis were subjected to SDS PAGE followed by western blot with CAF50452(ABC3A Antibody) at dilution of 1:1000

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.