

anti- DDIT3 antibody

Product Information

Catalog No.:	CAF50477
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information


Immunogen:	DNA-damage-inducible transcript 3
Synonyms:	CEBPZ, CHOP, CHOP 10, DDIT 3, DDIT3, GADD153
Observed MW:	19 kDa, 30 kDa
Uniprot ID :	P35638

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



Immunohistochemistry of paraffin-embedded human cervical cancer using CAF50477(CHOP antibody) at dilution of 1:50



K562 cells were subjected to SDS PAGE followed by western blot with CAF50477(DDIT3 antibody) at dilution of 1:1000

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.