

anti- RORC antibody

Product Information

Catalog No.:	CAF50478
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	RAR-related orphan receptor C
Synonyms:	NR1F3, Nuclear receptor ROR gamma, Nuclear receptor RZR gamma, RAR related orphan receptor C, RORC, RORG, RZR GAMMA, RZRG, TOR
Observed MW:	58 kDa
Uniprot ID :	P51449

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.



HT29 cells were subjected to SDS PAGE followed by western blot with CAF50478(RORC antibody) at dilution of 1:1000

ONLY

SAMPLE

This Antibody is for Research Use Only. Not for Diagnostic Procedures.

This is a sample Antibody manual only. Always refer to the hard copy manual included in the Antibody for your experiment.