

CAU21886 Mouse

Polyclonal Antibody to Klotho Beta (KLb)

Organism Species: Mus musculus (Mouse)

Instruction manual

FOR RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)



[PROPERTIES]

Source: Polyclonal antibody preparation

Host: Rabbit

Purification: Antigen-specific affinity chromatography followed by Protein A affinity

chromatography

Traits: Liquid

Concentration: 0.91mg/ml

UOM: 10µg(11µl)

Applications: WB; IHC; ICC; IP.

[IMMUNOGEN]

Immunogen: Recombinant KLb (Ser380~Phe506 (Accession # Q99N32)) expressed in E.coli

Accession No.: RPH756Mu02

[APPLICATIONS]

Western blotting: 0.5-2µg/mL;1:450-1800

Immunohistochemistry: 5-20µg/mL;1:45-180

Immunocytochemistry: 5-20µg/mL;1:45-180

Optimal working dilutions must be determined by end user.

[FORMULATION]

Form & Buffer: Supplied as solution form in 0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 4°C for frequent use.

Aliquot and store at -20°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the



expiration date under appropriate storage condition.

[IDENTIFICATION]






Figure. Western Blot; Sample:

Recombinant Klb, Mouse.



DAB staining on IHC-P;

Samples: Mouse Stomach Tissue



DAB staining on IHC-P;

Samples: Mouse Testis Tissue;

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.