

 Biomatik

 Tel:
 (519)
 489-7195,
 (800)

 836-8089
 Fax: (519)
 231-0140,
 (877)
 221-3515

 Email:
 info@biomatik.com
 http://www.biomatik.com

CAU25182

Polyclonal Antibody to Matrix Metalloproteinase 21 (MMP21) Organism Species: Homo sapiens (Human)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

[PRODUCT INFORMATION]

Immunogen: MMP21, Human Clonality: Polyclonal Host: Rabbit Immunoglobulin Type: IgG Purification: Affinity Chromatography. Applications: WB, ICC, IHC-P, IHC-F, ELISA Concentration: 200µg/mL UOM: 100µg

[IMMUNOGEN INFORMATION]

Immunogen: Recombinant MMP21 (Thr335~Phe547) with two N-terminal Tags, His-tag and T7-tag expressed in *E.coli*. **Accession No.:** RPB857Hu01

[ANTIBODY SPECIFITY]

The antibody is a rabbit polyclonal antibody raised against MMP21. It has been selected for its ability to recognize MMP21 in immunohistochemical staining and western blotting.



 Biomatik

 Tel:
 (519)
 489-7195,
 (800)

 836-8089
 Fax: (519)
 231-0140,
 (877)
 221-3515

 Email:
 info@biomatik.com
 http://www.biomatik.com

[APPLICATIONS]

Western blotting: 1:50-400 Immunocytochemistry in formalin fixed cells: 1:50-500 Immunohistochemistry in formalin fixed frozen section: 1:50-500 Immunohistochemistry in paraffin section: 1:10-100 Enzyme-linked Immunosorbent Assay: 1:100-200 Optimal working dilutions must be determined by end user.

[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH7.4, containing 0.02% NaN₃, 50% glycerol.

[QUALITY CONTROL]

Content: The quality control contains recombinant MMP21 (Thr335~Phe547) disposed in loading buffer.

Usage: 10uL per well when 3,3'-Diaminobenzidine(DAB) as the substrate.
 5uL per well when used in enhanced chemilumescent (ECL).
 Note: The quality control is specifically manufactured as the positive control.
 Not used for other purposes.

Loading Buffer: 100mM Tris(pH8.8), 2% SDS, 200mM NaCl, 50% glycerol, BPB 0.01%, NaN₃ 0.02%.

[STORAGE]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.


#1 Source for BioReagents, ELISA Kits and Custom Gene, Peptide, Protein & Antibody Services




Biomatik

Tel: (519) 489-7195, (800) 836-8089 Fax: (519) 231-0140, (877) 221-3515 Email: info@biomatik.com http://www.biomatik.com

[<u>IMAGES</u>]



Used in Western Blot, Sample: Recombinant MMP21, Human



Used in DAB staining on fromalin fixed paraffin- embedded liver tissue

#1 Source for BioReagents, ELISA Kits and Custom Gene, Peptide, Protein & Antibody Services