

SARS-CoV-2-S1 protein

Catalog No.: RPD50001

Application: SARS-CoV-2 kit

Source: Recombinant SARS-CoV-2-S1 protein obtained using mammalian cell

Structure: Consist of 695a.a. protein with C-terminal 6xHis tag.

Molecular weight: Glycosylated protein with molecular weight of ~118 kD

Purity: >95% as analyzed by SDS-PAGE

Buffer: 1 × PBS, pH7.4

Storage: Upon delivery aliquot and store at -80°C for up to 3 years. Avoid freeze/thaw cycles.