

1mg **RPU52716** 

**Recombinant Sestrin 1 (SESN1)** 

Organism Species: Mus musculus (Mouse)

Instruction manual

FOR RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)



### [PROPERTIES]

**Source:** Prokaryotic expression

Host: E.coli

Residues: Gly45~Thr492 Tags:

N-terminal His Tag Subcellular

Location: Nucleus Purity: >

85%

**Traits:** Freeze-dried powder

**Buffer formulation:** 20mM Tris, 150mM NaCl, pH8.0, containing 0.01% SKL, 5% Trehalose.

Original Concentration: 250µg/mL

**Applications:** Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 5.7

Predicted Molecular Mass: 55.7kDa

Accurate Molecular Mass: 56kDa as determined by SDS-PAGE reducing conditions.

#### [USAGE]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

#### [STORAGE AND STABILITY]

**Storage:** Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.


**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

#### [SEQUENCE]



				GHGPSR
FIPEKEMLQV	GSEDAQMHAL	FADSFAALGR	LDNITLVMVF	HPQYLESFLK
TQHYLLQMDG	PLPLHYRHYI	GIMAAARHQC	SYLVNLHVSD	FLHVGGDPKW
LNGLENAPQK	LQNLGELNKV	LAHRPWLITK	EHIEGLLKAE	EHSWSLAELV
HAVVLLTHYH	SLASFTFGCG	ISPEIHCDGG	HTFRPPSVSN	YCICDITNGN
HSVDEMQVNS	AGNASVSDSF	FEVEALMEKM	RQLQECREEE	EASQEEMASR
FEMEKRESMF	VFSSDDDEVT	PARDVSRHFE	DTSYGYKDFS	RHGMHVPTFR
VQDYCWEDHG	YSLVNRLYPD	VGQLIDEKFH	IAYNLTYNTM	AMHKDVDTSM
LRRAIWNYIH	CMFGIRYDDY	DYGEINQLLD	RSFKVYIKTV	VCTPEKVTKR
MYDSFWRQFK	HSEKVHVNLL	LIEARMQAEL	LYALRAITRY	MT

## [ IDENTIFICATION ]



# [ IMPORTANT NOTE ]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.